

2023 / 24

DENTAL HANDPIECE & ASSOCIATED ITEMS

Dentist's Office Reference



**Nakamura
Dental Handpiece
Mfg. Co., Ltd.**



E-TYPE CONNECTION HANDPIECE SET



ME-20MS

includes

- ES-30A** Straight Handpiece
- EC-20L** Latch Head Contra Angle
- ME-20M** 20,000rpm Airmotor



ES-30A

Straight Handpiece
1 year warranty



EC-20L

Contra Angle Handpiece
with Standard Latch Head
Warranty / Head : 3 months, Sheath : 6 months

※ *Contra Angle in the set can be exchanged with other models.
Please see page 5 for more information of other contra angles.*



ME-20M (4 holes)

ME-20B (2 holes)

20,000rpm Airmotor
1 year warranty



ME-80M (4 holes)

ME-80B (2 holes)

5,000rpm Airmotor
1 year warranty



STAR TITAN TYPE CONNECTION HANDPIECE SET



MS-10MHS

includes

- STS-30H** Straight Handpiece
- STCH-20L** Latch Head Contra Angle
- MS-10MH** 20,000rpm Airmotor



STS-30H

Straight Handpiece
1 year warranty



STCH-20L

Contra Angle Handpiece with Standard Latch Head
Warranty / Head : 3 months, Sheath : 6 months



STUNH-20L

Contra Angle Handpiece with Standard Latch Head
Warranty / Head : 3 months, Sheath : 6 months

※ *Contra Angle in the set can be exchanged with other models.
Please see page 6,7 for more information of other contra angles.*



MS-10MH

20,000rpm Airmotor
1 year warranty



MS-55MH

5,000rpm Airmotor
1 year warranty



MS-ADPH

4:1 Speed Reduction
Torque Multiplier
6 month warranty



AIRMOTOR HANDPIECE



MM-50M (5,000rpm)
MM-20M (20,000rpm)

Slim design, light weight, smooth run, minimum heat generation achieve less fatigue for its user and guarantee trouble free operation.



MD-50M (5,000rpm)
MD-20M (20,000rpm)

1 Year Warranty

U-TYPE CONTRA ANGLE FOR MM-50M / MM-20M / MD-50M / MD-20M



DC-10L
15,000rpm for CA Bur



DC-10FG
15,000rpm for FG Bur



DU-20L
20,000rpm for CA Bur
Spare Head : EH-20L



DU-30BL
30,000rpm for CA Bur
Spare Head : EH-30BL



DU-30TL
Quarter Turn For Endodontic
Spare Head : EH-30TL



DU-20FG
15,000rpm for FG Bur
Spare Head : EH-20FG



DU-35BL
30,000rpm for CA Bur
Spare Head : EH-35BL

3 Month Warranty




AIRMOTOR HANDPIECE FOR HYGIENIST

135°C
)))



MP-40M



MPS-40M

- The compact size with light weight is ergonomically designed to reduce fatigue for long time use.
- Can accept almost all disposable prophy angles
- 5,000rpm


1 Year Warranty

PROPHY ANGLE FOR U-TYPE NOSE

135°C
)))



RA-10PS
Regular Screw-in Propphy Angle



RA-10PK
Regular Snap-on Propphy Angle



RA-25PS
Sealed-type Screw-in Propphy Angle



RA-25PK
Sealed-type Snap-on Propphy Angle

3 Month Warranty



CONTRA ANGLE HANDPIECE



E-type



EC-20L
Standard Latch Head




EC-20ML
Miniature Latch Head



EC-20FG
Friction Grip Head



EC-30BL
Push Up, Ball Bearing Head



EC-35BL
Side Swing Latch, Ball Bearing Head



EC-30BLP
Push Button, Ball Bearing Head



CGM-20L
Standard Latch Head
4:1 Speed Reduction



CGM-20ML
Miniature Latch Head
4:1 Speed Reduction



CGM-20FG
Friction Grip Head
4:1 Speed Reduction



CGM-30TL
Quarter Turn Head for Endodontic
4:1 Speed Reduction



CGM-30BL
Push Up Latch, Ball Bearing Head
4:1 Speed Reduction




CGM-35BL
Side Swing Latch, Ball Bearing Head
4:1 Speed Reduction



CONTRA ANGLE HANDPIECE




CGM-20PS
Regular Prohphy, Screw-in
4:1 Speed Reduction



CGM-50PS
Sealed Prohphy Screw-in,
4:1 Speed Reduction




CGM-20PK
Regular Prohphy, Snap-on
4:1 Speed Reduction



CGM-50PK
Sealed Prohphy, Snap-on
4:1 Speed Reduction



CRM-20L
Standard Latch Head
16:1 Speed Reduction



CRM-30TL
Quarter Turn Head for Endodontic
16:1 Speed Reduction

HANDPIECE MAINTENANCE OIL SPRAY


YOSHIDA
Lubrication Oil Spray



Star Titan - type



STUNH-20L
Standard Latch Head



STUNH-20ML
Miniature Latch Head



CONTRA ANGLE HANDPIECE



Star Titan-type



STUNH-20FG
Friction Grip Head



STUNH-30BL
Push Up Latch, Ball Bearing Head




STUNH-30BLP
Push Button, Ball Bearing Head




STUNH-20PS
Regular Prohphy, Screw-in




STUNH-20PK
Regular Prohphy, Snap-on



STCH-20L
Standard Latch Head



STCH-20PS
Regular Prohphy, Screw-in



STUNH-30TL
Quarter Turn Head for Endodontic



STUNH-35BL
Side Swing Latch, Ball Bearing Head



STUNH-50PS
Sealed Prohphy, Screw-in



STUNH-50PK
Sealed Prohphy, Snap-on



STCH-30BL
Push Up Latch, Ball Bearing Head




CONTRA ANGLE HANDPIECE



Midwest-type




MWC-20L
Standard Latch Head



MWC-20ML
Miniature Latch Head



MWC-20FG
Friction Grip Head



MWC-30TL
Quarter Turn Head for Endodontic



MWC-30BL
Push Up Latch, Ball Bearing Head



MWC-35BL
Side Swing Latch, Ball Bearing Head



MWC-30BLP
Push Button, Ball Bearing Head




MWC-20PS
Regular Prohphy, Screw-in



MWC-50PS
Sealed Prohphy, Screw-in



MWC-20PK
Regular Prohphy, Snap-on



MWC-50PK
Sealed Prohphy, Snap-on



(👉) NEW PRODUCTS!

PUSH-BUTTON BALL-BEARING FOR FRICTION GRIP BUR




HEAD ATTACHMENT



EH-30FGP (14-teeth)
EHN-30FGP (12-teeth)

3 Month Warranty


CONTRA ANGLE HANDPIECE



E-Type
EC-30FGP (14-teeth)
ECN-30FGP (12-teeth)



Star-Titan Type
STUNH-30FGP



Midwest Type
MWC-30FGP

Head: 3 Month Warranty
Sheath: 6 Month Warranty



STRIPPING CONTRA ANGLE HANDPIECE



EG-50VT

E-type Connection
 1.3 mm Stroke Length
 4:1 Speed reduction for excellent torque
 Vertical reciprocating move



EG-10VT

E-type Connection
 0.4 mm Stroke Length
 4:1 Speed reduction for excellent torque
 Vertical reciprocating move



Vertical Reciprocating Head

EH-50VT

1.3mm Stroke Length

EH-10VT

0.4mm Stroke Length



AUTO STRIPPING BUR - Double Sided



Double Sided Strip is a thin strip with diamond coating on two sides.



PRODUCT CODE	DENSITY OF DIAMOND COATING ON THE BLADE (Unit : μm)	THICKNESS OF BLADE (Unit : mm)
001-015H	+ 15 (FINE)	+ 0.08
001-025H	+ 25	+ 0.12
001-040H	+ 40	+ 0.15
001-060H	+ 60	+ 0.19
001-090H	+ 90 (COARSE)	+ 0.29

Thickness of Blade is not same size as a gap between the teeth to strip.
 1μm + 1/1000 mm.



E-TYPE / MIDWEST-TYPE HANDPIECE SYSTEM



All the head attachments are 3 month warranty.



EH-20L
Latch-type Standard
for CA bur, 20,000rpm



EH-20ML
Latch-type Miniature
for CA bur, 20,000rpm



EH-20FG
For FG bur, Tool Type
20,000rpm



EH-30TL
Quarter Turn Move for
Endodontic for CA bur



EH-30BLP
Push Button, Ball-bearing
for CA bur, 30,000rpm



EH-30BL
Push-up Latch,
Ball-bearing
for CA bur, 30,000rpm



EH-35BL
Latch-type,
Ball-bearing
for CA bur, 30,000rpm




EH-20PS
Regular Screw-in
Prophy



EH-20PK
Regular Snap-on
Prophy



EH-50PS
Sealed-type, Screw-in Prophy



EH-50PK
Sealed-type, Snap-on Prophy



SHS-EC
E-type Contra Angle Sheath
6 month warranty



SHM-E04
4:1 Speed Reduction
E-type Contra Angle Sheath
6 month warranty



SHM-E16
16:1 Speed Reduction
E-type Contra Angle Sheath
6 month warranty



SHS-MWC
Midwest-type
Contra Angle Sheath
6 month warranty



E-TYPE / STAR TITAN TYPE HANDPIECE SYSTEM



All the head attachments are 3 month warranty.



MHN-20L
For CA bur, 20,000rpm



MHN-30BL
Ball-Bearing for CA bur
30,000rpm



MHN-20PS
Regular prophy
Screw-on



SHS-STCH
Contra Angle Sheath
6 month warranty




EHN-20L
Latch-type Standard
for CA bur, 20,000rpm



EHN-20ML
Latch-type, Miniature
for CA bur, 20,000rpm



EHN-20FG
For FG bur, Tool Type
20,000rpm



EHN-30TL
Quarter Turn Move for
Endodontic for CA bur



SHS-STUNH
Contra Angle Sheath
6 month warranty




SHS-ECN
E-type Contra Angle Sheath
6 month warranty



EHN-30BLP
Push Button, Ball-bearing
for CA bur, 30,000rpm



EHN-30BL
Push-up Latch, Ball-bearing
for CA bur, 30,000rpm



EHN-35BL
Latch-type, Ball-bearing
for CA bur, 30,000rpm



SHN-E04
4:1 Speed Reduction
E-type Contra Angle Sheath
6 month warranty




EHN-20PS
Regular Screw-in Prophy



EHN-20PK
Regular Snap-on Prophy



EHN-50PS
Sealed-type, Screw-in Prophy



EHN-50PK
Sealed-type, Snap-on Prophy



SHN-E16
16:1 Speed Reduction
E-type Contra Angle Sheath
6 month warranty



STRAIGHT HANDPIECE



ES-30A
E-Type



ESG-30AR
E-Type 4:1 Speed Reduction



STS-30H
Star-Titan type



MWS-30A
Midwest-type

1 Year Warranty

AIRTURBINE HANDPIECE



Push-Button Chuck Type



TCP-M01
Standard Size Head, 400,000rpm



CRT-M01
Cartridge for TCP-M01
3 month warranty

Manual Chuck Type



TCH-M02
Standard Size Head, 400,000rpm



CRT-M02
Cartridge for TCH-M02
3 month warranty



LABORATORY HANDPIECE



TS-300ZT
Regular Airturbine Handpiece
Warranty / Cartridge : 3 months, Body : 1 year



CRT-300Z
Cartridge for TS-300ZT
3 month warranty



TS-100LT
Lube-Free Airturbine Handpiece
Warranty / Cartridge : 3 months, Body : 1 year



CRT-100L
Cartridge for TS-100LT
3 month warranty



LABORATORY MICROMOTOR SET

LMS-20S

includes

LBM-100 2,000 - 20,000rpm E-type Micromotor

CNT-001 Controller

ES-30A E-type Straight Handpiece

AC Adaptor

1 year warranty

※ Battery pack is not included in the set.

