

FINENESS

Ultra Thin, Light Weight Resin Forming Instrument "FINENESS" Lines

The instrument coated by titanium nitride, is comfortable to use. Ultra thin tip achieves easy access to fill and form resin in area that has traditionally been difficult to reach.

Ultra Thin
0.17mm

Light
Weight
15g

Hollow stainless steel handle is adopted. You have great selection among all the lines of easy-to-use instruments.

Ultra Thin, Light Weight 4 Different Fineness

Fineness 1




Fineness 2



Fineness 3



Fineness MI



Fineness

Please see the back side for product specification of Resin Forming Instrument "FINENESS" Lines.

Resin forming is made easy in area that has traditionally been difficult to reach.

Resin Forming Instrument with Great Comfort & Convenience To Use.

- Adopted titanium nitride coating prevents resin from sticking on surface.
- Hollow stainless steel handle with light weight (15g) reduces your hand fatigue.
- Ultra thin tip enables you to fill and form resin in area that has traditionally been difficult to reach.

Ideal for filling and forming composite / flowable resin.



Flat Face




Anterior Tooth Side



Tooth Bottom Area

6 points of Fineness Lines

Check



1.5 as wide as conventional instrument.



2. Wide & short type is available, too.




3. Resin filling for triangular hollow.



4. Narrow area between teeth can be accessed.



5. For tiny filling, removal of excessive material, use with microscope.



6. Bent form with great flexibility to securely touch tooth face.

Please contact to
HEAD DENTAL CORPORATION
17972 Sky Park Circle Suite J
Irvine, CA 92614-6409
Phone 949-474-0176 / FAX 949-474-1736
charlie@headdental.com

Manufactured by
SAITO DENTAL MFG. Co.,LTD
2-56-8 Horikiri Katsushika-ku TOKYO JAPAN
Phone (81-3)3691-2219 / FAX (81-3)3691-9487
hiro@saito-dental.co.jp